EdgeOne-R Utopia R-Glide Instructions For Use

Ref.: 30006123-A

Character	N	L	Stop	Req.	RPM	Angles	Maximum recommended torque
R-Glide	14	50	21-30 Hz	White	100-400 rpm	CE marked	4.0

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- 5. Clinical claims**
EdgeOne-R Utopia R-Glide files, in normal conditions. Clinical performance:
- Shaping ability: adapted debris removal and cutting efficiency
 - Respect of canal anatomy
 - Geometry adapted to irrigation
 - Maintains the apical foramen in its original position
- Usability:
- Secure progression and limit the interference
- Safety:
- Limited risk of breakage
 - Limited risk of shaping aberration
 - Limited screwing effects

- 6. Characteristics**
Instrument in reciprocating motion in type 4 under EN ISO 3630-1: 2019 standard: rotation angles are explained in the column «CW» (Counter-clockwise side) and column «CW» (Clockwise side).

- Material of operative part: Nitinol.
Use in connection with an endodontic contra-angle according to EN ISO 1797: 2017 (Type 1).
Number of uses: single use.
Medical device class according to Directive 93/42/EEC and MDR 2017/745: Ia.
7. Protocol
General information on the protocol for use:
After opening the access cavity, sodium hypochlorite must be constantly used and renewed before and after each passage of instruments. The sodium hypochlorite must be placed using the syringe in order to be injected as close as possible to the apical region.
Operating procedure:
EdgeOne-R Utopia R-Glide is introduced into the canal with a slow and regular movement without any apical pressure until the working length is reached.
If an ISO size #10 hand file will go to working length, then make a mechanized glidepath with EdgeOne-R Utopia R-Glide up to the working length and insert again the ISO size #10 hand file up to the working length.
If an ISO size #10 hand file will not go to working length, then initially make a mechanized glidepath with EdgeOne-R Utopia R-Glide up to the penetration length of the ISO size #10 hand file. Then electronically determine the working length with an ISO size #10 hand file. And then make a mechanized glidepath with EdgeOne-R Utopia R-Glide up to the working length and insert again the ISO size #10 hand file up to the working length.

- 8. Reprocessing instructions**
Not applicable.

- 9. Storage and transport conditions**
After use, instruments must be placed in a secure container, used to collect cutting or sticking instruments (like needles or disposable bistouries) as per good dentistry practices.

- 10. Disposal**
After use, instruments must be placed in a secure container, used to collect cutting or sticking instruments (like needles or disposable bistouries) as per good dentistry practices.

11. Symbols			
	NITI raw material		Root canal preparation
	Quantity		Assortment
	"CLASSICS" handle		Do not use if packaging is damaged
	Medical device		

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- 8. Instructions de reprocessing**
Not applicable.

- 9. Instruction de stockage et transport**
Non applicable.

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