INSTRUCTIONS FOR USE – EDGE X7 UTOPIA

EN	INSTRUCTIONS FOR USE	2
AR	تعليمات الاستخدام	5
BG	ИНСТРУКЦИИ ЗА УПОТРЕБА	8
CS	NÁVOD K POUŽITÍ	11
DA	BRUGSANVISNING	14
DE	GEBRAUCHSANWEISUNG	17
EL	ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ	20
ES	INSTRUCCIONES DE USO	24
ET	KASUTUSJUHEND.....	27
FI	KÄYTTÖOHJEET	30
FR	NOTICE D'UTILISATION.....	33
HR	UPUTE ZA UPORABU	36
HU	HASZNALATI UTASITAS.....	39
IT	ISTRUZIONI PER L'USO	42
LT	NAUDOJIMO INSTRUKCIJA	45
LV	LIETOŠANAS INSTRUKCIJAS	48
NL	GEBRUIKSAANWIJZING	51
NO	BRUKSANVISNING	54
PL	INSTRUKCJA UŻYWANIA.....	57
PT	INSTRUÇÕES DE UTILIZAÇÃO	60
RO	INSTRUCȚIUNI DE UTILIZARE	63
SK	NÁVOD NA POUŽITIE.....	66
SL	NAVODILA ZA UPORABO	69
SV	BRUKSANVISNING	72
TR	KULLANIM TALIMATLARI	75

INSTRUCTIONS FOR USE


EN

FOR DENTAL USE ONLY


1) GENERAL DESCRIPTION

Composition:

The instrument is made of nickel-titanium alloy and the endo stop is made of silicone.

Tip size and taper are identified by two color markings: the wide colored ring on the handle indicates the tip diameter, while the stopper and the narrow colored ring indicate the taper.

Parameters for use:

These instruments are compatible with endodontic motors compliant with ISO 1797.
Recommended speed: 500 rpm continuous rotation; torque limit: 410 g·cm (4 N·cm).

2) INTENDED USE

The Edge X7 Utopia™ endodontic instruments are intended for root canal treatment.

3) INDICATIONS FOR USE

Edge X7 Utopia™ are indicated in case of pulp or periapical diseases that require root canal treatment.
All instruments are designed for preparation and shaping of root canals.

4) CONTRAINDICATIONS AND ADVERSE REACTIONS

None known.

5) WARNINGS AND PRECAUTIONS

- These instruments have not been tested on children, pregnant nor breastfeeding women.
- These endodontic instruments must not be used on individuals with a known allergic sensitivity to nickel-titanium and silicone.
- Single use instruments must not be reprocessed for subsequent use nor be reused. This could compromise their integrity.
- Each instrument may be used to treat up to four (4) canals during a single procedure. Discard after use.
- For your own safety, wear personal protective equipment (masks, gloves, and goggles) during the treatment and reprocessing of the instrument.
- Do not use if sterile packaging seal is broken, damaged or wet.
- Do not use the instruments after expiration date.
- Reprocess the instrument appropriately before disposal.
- Carefully read the label, marking on the packaging and this instruction for use to ensure the correct identification and use of device.
- Use a dental dam when using endodontic instruments to avoid inhalation or ingestion by the patient.
- Take several radiographs from different angles to determine the anatomy of the root canals (length, width and curvature).
- Always inspect the instrument(s) before use, and discard (it) them if there are any visible defect(s).
- Cracks, deformations, signs of corrosion, loss in color or marking are signs that the device is no longer able to achieve the required performance level and should be discarded.
- Regularly clean flutes and inspect instrument during use, and discard if it shows any signs of distortion or wear.
- If the packaging has been opened, damaged, become wet or the expiry date has passed, the sterile state of the instruments inside the packaging is not guaranteed. Discard the device(s) that are single use and consult the reprocessing instructions available on EdgeEndo's website for multiple use device(s).
- All serious events occurring in relation to the product must be reported to the manufacturer and the competent authority according to local regulations.

6) PERFORMANCE CHARACTERISTICS

Edge X7 Utopia™ rotary files are compliant with all applicable testing requirements of the ISO 3630 series of standards.

7) PROTOCOL FOR USE

1) Simplified Technique Guide

- a. Create straight line access to canal orifices.
- b. Locate canals and explore using stainless steel hand instruments.
- c. Irrigate before each hand or rotary file with liquid NaOCl or 17% EDTA.
- d. Using a #10 K file or larger file determine the working length with an apex locator and/or radiographs in each canal. This will also form a glide path in each canal.
- e. Use #17/04 file in one or more passes, alternating with small-sized hand files if necessary, until working length is reached. If more coronal flare is desired it can be achieved by incorporating the EdgeGlidePath file used in a brushing motion.
- f. Next use #25/04 to working length passively; if instrument has not reached working length use smaller additional shaping instrument #20/04 to working length alternating between #20/04 and #25/04.

INSTRUCTIONS FOR USE

EN

- g. After #25/04 reaches working length if clinician desires a larger apical or canal shape, take a #25/06 or other additional larger 04 or 06 instruments (#30, #35, #40, etc...).

2) Disinfecting

- After each canal is fully shaped, rinse the canals for 1 minute with 17% Liquid EDTA to remove the canal Smear Layer.
- Rinse the canals for 5 minutes with 5% NaOCl to remove debris and bacteria.
- Rinse the canals for 1 minute with 17% Liquid EDTA to rinse out the 5% NaOCl.
- Rinse the canals for 5 minutes with 2% chlorohexidine or EDTA to kill bacteria.

3) Obturation

- When using a thermal carrier system, use size verifiers to determine the proper sized carrier.
- When using a master gutta percha cone that matches the largest file taken to length, remember sometimes you may need to drop down in cone tip size if the corresponding gutta percha to your final rotary file does not go to length.

8) CONDITIONS AND TERM OF STORAGE, SHELF-LIFE

The shelf-life is 5 years for sterile instruments. The use-by date is indicated on the packaging. Sterility is preserved and guaranteed for the whole shelf life unless the packaging is broken, or the storage conditions are not respected.

9) DISPOSAL

When a device reaches the end of its life, make sure that it is discarded in accordance with the applicable laws and regulations.

10) SYMBOLS

For explanation of symbols for IFUs and labels, please consult the document “Glossary of the symbols used by EdgeEndo” available on EdgeEndo’s website.

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Henry Schein Services GmbH,