INSTRUCTIONS FOR USE – EDGEONE FIRE

EN	INSTRUCTIONS FOR USE	2
AR	تعليمات الاستخدام	6
BG	ИНСТРУКЦИИ ЗА УПОТРЕБА.....	9
CS	NÁVOD K POUŽITÍ	13
DA	BRUGSANVISNING.....	17
DE	GEBRAUCHSANWEISUNG	21
EL	ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ	25
ES	INSTRUCCIONES DE USO.....	29
ET	KASUTUSJUHEND	33
FI	KÄYTTÖOHJE	37
FR	NOTICE D’UTILISATION	41
HR	UPUTE ZA UPORABU	45
HU	HASZNÁLATI UTASÍTÁS	49
IT	ISTRUZIONI PER L’USO.....	53
LT	NAUDOJIMO INSTRUKCIJA.....	57
LV	LIETOŠANAS INSTRUKCIJA	61
NL	GEBRUIKSAANWIJZING.....	65
NO	BRUKSANVISNING	69
PL	INSTRUKCJA UŻYWANIA	73
PT	INSTRUÇÕES DE UTILIZAÇÃO	77
RO	INSTRUCȚIUNI DE UTILIZARE.....	81
SK	NÁVOD NA POUŽITIE	85
SL	NAVODILA ZA UPORABO	89
SV	BRUKSANVISNING	93
TR	KULLANIM TALIMATLARI	97

INSTRUCTIONS FOR USE

EN


FOR DENTAL USE ONLY

INSTRUCTIONS FOR USE – EDGEONE FIRE

1) GENERAL DESCRIPTION

Composition:

The instrument is made of nickel-titanium alloy and the endo stop is made of silicone.

The tip size is indicated by the color band. The taper is variable.

These instruments are compatible with endodontic motors compliant with ISO 1797.

Recommended speed: 350 rpm Reciprocating motion; torque limit: 408–530 g·cm (4.0–5.2 N·cm).

Name	Size	Color
Small	020	Yellow
Primary	025	Red
Medium	035	Green
Large	045	White

2) INTENDED USE

The EdgeOne Fire™ endodontic instruments are intended for root canal treatment.

3) INDICATIONS FOR USE

EdgeOne Fire™ are indicated in case of pulp or periapical diseases that require root canal treatment. All instruments are designed for preparation and shaping of root canals.

4) CONTRAINDICATIONS AND ADVERSE REACTIONS

None known.

5) WARNINGS AND PRECAUTIONS

- These instruments have not been tested on children, pregnant nor breastfeeding women.
- These endodontic instruments must not be used on individuals with a known allergic sensitivity to nickel-titanium and silicone.
- Single use instruments must not be reprocessed for subsequent use nor be reused. This could compromise their integrity.
- Each instrument may be used to treat up to four (4) canals during a single procedure. Discard after use.
- For your own safety, wear personal protective equipment (masks, gloves, and goggles) during the treatment and reprocessing of the instrument.
- Do not use if sterile packaging seal is broken, damaged or wet.
- Do not use the instruments after expiration date.
- Reprocess the instrument appropriately before disposal.
- Carefully read the label, marking on the packaging and this instruction for use to ensure the correct identification and use of device.
- Use a dental dam when using endodontic instruments to avoid inhalation or ingestion by the patient.
- Take several radiographs from different angles to determine the anatomy of the root canals (length, width and curvature).
- Always inspect the instrument(s) before use, and discard (it) them if there are any visible defect(s).
- Cracks, deformations, signs of corrosion, loss in color or marking are signs that the device is no longer able to achieve the required performance level and should be discarded.
- Regularly clean flutes and inspect instrument during use, and discard if it shows any signs of distortion or wear.
- The packaging of the instruments should be carefully examined before opening (packaging integrity, no humidity, no visual degradation of instruments) to ensure that the packaging's integrity has not been compromised during transport and storage.
- If the packaging has been opened, damaged, become wet or the expiry date has passed, the sterile state of the instruments inside the packaging is not guaranteed. Discard the device(s) that are single use and consult the reprocessing instructions available on EdgeEndo's website for multiple use device(s).
- All serious events occurring in relation to the product must be reported to the manufacturer and the competent authority according to local regulations.

6) PERFORMANCE CHARACTERISTICS

EdgeOne Fire™ rotary files are compliant with all applicable testing requirements of the ISO 3630 series of standards.

7) PROTOCOL FOR USE

1) Size Selection:

- a. If the #10 hand file was tight, use the EdgeOne Fire™ Small.
- b. If the #10 hand file was easy but the #15 hand file was tight, use the EdgeOne Fire™ Medium.
- c. If both the #10 and #15 hand files were easy, use the EdgeOne Fire™ Large.

2) Glide Path:

- Create a glide path and determine the working length prior to EdgeOne Fire™ file use by negotiating all root canals to their terminus with stainless steel #10 and #15 hand files and a lubricant.

3) Canal Shaping and Cleaning:

- a. Place the selected EdgeOne Fire™ file into the handpiece and set it for reciprocation motion.
- b. With lubricant in the canal and light apical pressure, use a gentle inward pecking motion, advancing the file 2 to 3 mm, then lifting up 1 to 2 mm. Keep repeating this motion to passively advance the EdgeOne Fire™ file until it does not easily progress.
- c. Remove the EdgeOne Fire™ file from the canal, remove debris and inspect the file, irrigate and recapitulate with a #10 hand file and check patency.
- d. Repeat steps c. and d. until the EdgeOne Fire™ file is at the working length. If, after repeated attempts, the EdgeOne Fire™ file does not seem to be advancing any further, use the next smaller EdgeOne Fire™ file size and finish the canal.
- e. Apically gauge the size of the foramen with a hand file of the same tip size as the EdgeOne Fire™ file taken to length. If the gauging hand file is a snug fit, the preparation is finished. If it is loose, use the next EdgeOne Fire™ file to finish the preparation.

4) Final Irrigation and Drying:

Perform a final irrigation protocol according to clinical preference. Dry the canal with paper points that correspond to the final file size taken to working length.

5) Obturation:

Obturate the canal system using a matching gutta-percha cone or thermal carrier system.

8) CONDITIONS AND TERM OF STORAGE, SHELF-LIFE

The shelf-life is 5 years for sterile instruments. The use-by date is indicated on the packaging. Sterility is preserved and guaranteed for the whole shelf life unless the packaging is broken, or the storage conditions are not respected.

9) DISPOSAL

When a device reaches the end of its life, make sure that it is discarded in accordance with the applicable laws and regulations.

INSTRUCTIONS FOR USE

EN

10) SYMBOLS

For explanation of symbols for IFUs and labels, please consult the document “Glossary of the symbols used by EdgeEndo” available on EdgeEndo’s website.

 **Henry Schein Services GmbH,**
Siemensstraße 1,
63263 Neu-Isenburg, Germany
Importer for the European Union

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