INSTRUCTIONS FOR USE – EDGETAPER & EDGETAPER PLATINUM

EN	INSTRUCTIONS FOR USE	2
AR	تعليمات الاستخدام	6
BG	ИНСТРУКЦИИ ЗА УПОТРЕБА	10
CS	NÁVOD K POUŽITÍ	14
DA	BRUGSANVISNING	18
DE	GEBRAUCHSANWEISUNG	22
EL	ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ	26
ES	INSTRUCCIONES DE USO	30
ET	KASUTUSJUHEND	34
FI	KÄYTTÖOHJEET	38
FR	NOTICE D'UTILISATION	42
HR	UPUTE ZA UPORABU	46
HU	HASZNALATI UTASITAS	50
IT	ISTRUZIONI PER L'USO	54
LT	NAUDOJIMO INSTRUKCIJA	58
LV	LIETOŠANAS INSTRUKCIJAS	62
NL	GEBRUIKSAANWIJZING	66
NO	BRUKSANVISNING	70
PL	INSTRUKCJA UŻYWANIA	74
PT	INSTRUÇÕES DE UTILIZAÇÃO	78
RO	INSTRUCȚIUNI DE UTILIZARE	82
SK	NÁVOD NA POUŽITIE	86
SL	NAVODILA ZA UPORABO	90
SV	BRUKSANVISNING	94
TR	KULLANIM TALIMATLARI	98

INSTRUCTIONS FOR USE


EN

FOR DENTAL USE ONLY

1) GENERAL DESCRIPTION

Composition:

The instrument is made of nickel-titanium alloy and the endo stop is made of silicone.

Instruments are identified using both the colored ring on the handle and the stopper color, as specified below:

Name	Size	Color on handle	Stopper Color	Taper	Shank Rings
S1	017	Purple	Purple	Variable	1
SX	019	No color	Yellow	Variable	0
S2	020	White	White	Variable	1
F1	020	Yellow	Yellow	Variable	1
F2	025	Red	Red	Variable	1
F3	030	Blue	Blue	Variable	1
F4	040	Black	Black	Variable	2
F5	050	Yellow	Black	Variable	2

No rings

1 ring

2 rings

These instruments are compatible with endodontic motors compliant with ISO 1797.

Recommended speed: 300-500 rpm continuous rotation; torque limit: 300 g·cm (2.9 N·cm).

2) INTENDED USE

The EdgeTaper™ and EdgeTaper Platinum™ endodontic instruments are intended for root canal treatment.

3) INDICATIONS FOR USE

EdgeTaper™ and EdgeTaper Platinum™ are indicated in case of pulp or periapical diseases that require root canal treatment. All instruments are designed for preparation and shaping of root canals.

4) CONTRAINDICATIONS AND ADVERSE REACTIONS

None known.

5) WARNINGS AND PRECAUTIONS

- These instruments have not been tested on children, pregnant nor breastfeeding women.
- These endodontic instruments must not be used on individuals with a known allergic sensitivity to nickel-titanium and silicone.
- Single use instruments must not be reprocessed for subsequent use nor be reused. This could compromise their integrity.
- Each instrument may be used to treat up to four (4) canals during a single procedure. Discard after use.
- For your own safety, wear personal protective equipment (masks, gloves, and goggles) during the treatment and reprocessing of the instrument.
- Do not use if sterile packaging seal is broken, damaged or wet.
- Do not use the instruments after expiration date.
- Reprocess the instrument appropriately before disposal.
- Carefully read the label, marking on the packaging and this instruction for use to ensure the correct identification and use of device.
- Use a dental dam when using endodontic instruments to avoid inhalation or ingestion by the patient.
- Take several radiographs from different angles to determine the anatomy of the root canals (length, width and curvature).
- Always inspect the instrument(s) before use, and discard (it) them if there are any visible defect(s).
- Cracks, deformations, signs of corrosion, loss in color or marking are signs that the device is no longer able to achieve the required performance level and should be discarded.
- Regularly clean flutes and inspect instrument during use, and discard if it shows any signs of distortion or wear.
- If the packaging has been opened, damaged, become wet or the expiry date has passed, the sterile state of the instruments inside the packaging is not guaranteed. Discard the device(s) that are single use and consult the reprocessing instructions available on EdgeEndo's website for multiple use device(s).
- All serious events occurring in relation to the product must be reported to the manufacturer and the competent authority according to local regulations.

6) PERFORMANCE CHARACTERISTICS

EdgeTaper™ and EdgeTaper Platinum™ rotary files are compliant with all applicable testing requirements of the ISO 3630 series of standards.

7) PROTOCOL FOR USE

1) Access and Canal location:

Create straight-line access to canal orifices. Locate canals and explore using stainless steel hand instruments.

2) Irrigation:

Irrigate the canal with an appropriate irrigating solution (e.g. sodium hypochlorite or EDTA) before each hand or rotary instrument.

3) Working Length and Glidepath:

- Determine working length using an apex locator and radiographs.
- Establish a glidepath to working length with a stainless-steel ISO #10 file or larger.
- Confirm canal patency with a #10 file.
- A mechanical glidepath such as EdgeGlide path instrument may be used according to the manufacturer's instructions.

4) Canal Shaping Using EdgeTaper™ and EdgeTaper Platinum™:

Straight/Wider Canals

- Use S1 to working length.
- Irrigate and recapitulate with a #10 hand file to working length.
- Use S2 to working length.
- Irrigate and recapitulate with a #10 hand file to working length.
- Use F1 to working length.
- Irrigate and recapitulate with a #10 hand file to working length.
- If a larger file is needed, then use F2 and F3 to working length with irrigation.

Curved/Narrow/Long Canals:

- Take a #10 hand file 2/3 down to the canal.
- A #15 hand file may also be used 2/3 down the canal, if needed.
- Take SX to the length of the hand files.
- Take a #10 hand file to estimated working length.
- Determine working length using an apex locator and radiographs.
- Confirm canal patency with a #10 file.
- Use S1 to working length.
- Irrigate, recapitulate, and check patency with a #10 hand file.
- Use S2 to working length.
- Irrigate, recapitulate, and check patency with a #10 hand file.
- Use F1 to working length.
- Irrigate, recapitulate, and check patency with a #10 hand file.
- If a larger file is needed, then use F2 and F3 to working length with irrigation.

5) Final Irrigation and Drying:

Perform a final irrigation protocol according to clinical preference. Dry the canal with paper points that correspond to the final file size taken to working length.

6) Obturation:

Obturate the canal system using a matching gutta-percha cone or thermal carrier system.

INSTRUCTIONS FOR USE

EN

8) CONDITIONS AND TERM OF STORAGE, SHELF-LIFE

The shelf-life is 5 years for sterile instruments. The use-by date is indicated on the packaging. Sterility is preserved and guaranteed for the whole shelf life unless the packaging is broken, or the storage conditions are not respected.

9) DISPOSAL

When a device reaches the end of its life, make sure that it is discarded in accordance with the applicable laws and regulations.

10) SYMBOLS

For explanation of symbols for IFUs and labels, please consult the document “Glossary of the symbols used by EdgeEndo” available on EdgeEndo’s website.

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Henry Schein Services GmbH,